



Excellence in Human Service
and Technology

Echo Sounder ES-101EX 2MHz FHR Model

HAND-HELD FETAL DOPPLER WITH LCD

Detects fetal heartbeat and displays the rate.

ACOUSTIC FEEDBACK CANCELLER

For clear sounds with built-in speaker & transducer.

EASY & SIMPLE OPERATION

One-handed operation includes freezing heart rate for notation.

CLINICAL APPLICATIONS

- Detection of fetal heart sounds
- Determination and monitoring of fetal heart rate
- Verification of fetal heart viability by patient trauma

DISPLAY SAMPLES



SPECIFICATIONS

Power:	DC 3V, 2 AA alkaline batteries
Battery life:	5 hours or more
Automatic shut-off:	3 min. with no signal
Speaker output:	Max. 250 mW
Frequency:	2.25 MHz
Heart rate range:	60 - 240 bpm
LCD display:	Heart rate, sound volume, low battery indicator
Output:	Headset output
Dimensions:	45 (W) x 135 (D) x 30 (H) mm
Weight:	Approx. 150 grams, including batteries
Electrical safety:	Conforms to IEC-60601-1

*Specifications subject to change.



101EX2 CL E S-10
Printed in Japan

Authorised Importer & Distributor:
ZOHRA WAHAB TRADING COMPANY PTE LTD, SINGAPORE
URL: www.zowamedical.com.sg / Tel: (65)-6743 6050 / 6747 0329