

AUTO - INFLATE DIGITAL BLOOD PRESSURE MONITOR

MODEL: KD-202F



KEY SPECIAL FEATURES:

- *Talking :Pure American pronunciation broadcast the measurement results*
- *Using the large LCD screen (97*87mm)with LED backlight*
- *Using the TOUCH PAD KEY*
- *Automatic measurement of systolic, diastolic and pulse with time & date*
- *Intelligent inflate, Reduces the oppression on the arm*
- *80 measurement results with time & date*
- *Automatically turns off to save power.*
- *CE, ISO13485 certification*
- *Comes with Adapter*

SPECIFICATIONS

Measurement method:	Oscillometric system
Display type:	Digital LCD 97mm [L] x 87mm [W]
Measuring range:	Pressure 0-300 mmHg Pulse 30~200 p/min
Measuring accuracy:	Pressure $\leq\pm 3$ mmHg Pulse $\leq 5\%$
Operating environment:	Temperature 10-40 \square Relative humidity $\leq 80\%$
Power source:	4 x AA batteries or AC/DC adapter
Dimension:	150mm [L] x 110mm [W] x 65mm [H]
Weight:	About 370g
Cuff size:	520mm [L] x 135mm [W]
Cuff range: size:	220mm [L] - 320mm [W]

PACKING DETAILS:

Color box size:	128mm [L] x 104 mm [W] x163mm [H]
Carton size:	555mm [L] x 355 mm [W] x300 mm [H]
Packing:	20pcs/ctn