


DS-137 Manual-inflate Digital Blood Pressure Monitor

for upper arm measurement



some of our first digital blood pressure monitors *currently not for sale

DS-137 is a digital blood pressure monitor that holds the classic style of those for home use – manual inflation of the cuff. Manual-inflate digital blood pressure monitor was developed almost 20 years ago. Although the style of measurement is classic, DS-137 system carries the most advanced technology to achieve accurate measurement of blood pressure as well as the iconic indication assisting measuring procedure. That is to say, DS-132 is the product in which 20 years of knowledge and technology are concentrated.



Product specifications


Model	DS-137
Principle of blood pressure measurement	oscillometric method
Indicator	11 digits liquid crystal display
Pressure indication range	0 to 300 mmHg
Measuring range	50 to 250 mmHg (systolic) 40 to 180 mmHg (diastolic) 40 to 160 pulses/min (pulse rate)
Accuracy	± 3 mmHg (blood pressure) ± 5 % of reading (pulse rate)
Inflation	squeeze bulb
Deflation	constant air release valve
Exhaust	manual rapid exhaust valve
Power supply	4 pcs. 1.5 volt AA (R6P/LR6) batteries
Power consumption	0.06W (max.)
Memory	30 measurement results and their average
Operating environment	+10°C to +40°C, 85 % relative humidity or below
Storage environment	-5°C to +50°C, 85 % relative humidity or below
Applicable arm circumference	22.0 to 32.0 cm
Dimensions	163 (W) x 180 (D) x 83 (H) mm, placed on stand
Weight	approximately 180 g, without batteries
Accessories	instruction manual, 4 AAA batteries



simultaneous measurement result display



easy operation with indication marks



stand tilting unit for better view with cuff pocket - Not included / Optional

Specifications are subject to change without prior notice due to improvements in performance and quality.

NISSEI 日本精密測器株式会社
NIHON SEIMITSU SOKKI CO., LTD. 2508-13 Komochi Kitagunma Gunma Japan
EC-Representative:
JPI Trading GmbH, Holzbrücke 7, 20459 Hamburg, Germany ☎ 040-82 22 80 510



web site <http://www.nissei-jp.com>

Importer/Distributor: **ZOWA**
(Since 1977)

Zohra Wahab Trading Company Pte Ltd, Singapore
TEL: (65)-6743 6050 / 6747 0329
URL: www.zowamedical.com.sg