



Japan (Since 1950)

HT-1500 ONE HAND PALM ANEROID GAUGE SPHYGMOMANOMETER



- Palm-style integrated aneroid sphygmomanometer
- Suitable for left or right handed use
- Compact and lightweight
- One touch trigger deflation
- Supplied with standard adult cuff and storage case

SPECIFICATIONS

- Pressure indication range : 0~300 mmHg
- Accuracy of cuff pressure display : ± 3 mmHg
- Applicable arm circumference : approximately 23 to 32 cm
- Operating Environment : +15°C to +25°C, 85% relative humidity or below
- Storage Environment : -20°C to +70°C, 85% relative humidity or below
- Size : 67×158.5×86.3 mm (W, D, H)
- Weight : 143 g

Importer/Distributor:

Zohra Wahab Trading Co Pte Ltd
629 Aljunied Road, #08-02, Cititech Industrial Building, Singapore 389838
Tel: 65-67436050 / 67470329 / Fax: 67471806
Website: www.zowamedical.com.sg