

DIGITAL HEARING AID

READY-MADE IN THE EAR TYPE



• LEFT EAR



• RIGHT EAR

- MILD TO MODERATE TRIMMER-TYPE DIGITAL
- NOISE REDUCTION & WATER-RESISTANT COATING
- ZINC AIR BATTERY 312 / PR 41 (1.4V)

HC-A1 Mi Mo

● Ready-made in-the-ear digital hearing aid for people with hearing disability of up to "Moderate"

There is no need to measure the shape of ear for fitting. The HC-A1 can be used immediately.

● Trimmer-type digital hearing aid

There is no need to use fitting software for making adjustments. Adjustments can be made easily by turning the two trimmer dials (OPC and TONE H).

● Right-ear and left-ear units are available.

● Nonlinear amplification

The HC-A1 amplifies low-level sounds clearly without causing loud sounds to echo.

● Adaptive feedback canceller (AFBC)

The AFBC automatically suppresses disturbing high-pitched noise.

● Noise reduction (NR)

The noise reduction function performs multi-channel ambient sound analysis and reduces noise accordingly.

Importer & Distributor:

ZOWA

TEL: (65)-6743 6050 / 6747 0329