










FITTING GUIDE

Classification of loss	Mi Mild	Mo Moderate	Se Severe	Pr Profound
Hearing Level (dB)	25 30 35	40 50 60	70 80	90 100
Ability to Understand Speech	<p>Difficulty understanding soft-spoken speech</p> 	<p>Understanding of speech at 1 to 1.5 m. Difficulty in group and classroom discussion</p> 	<p>Loud speech may be understood at 30cm from ear</p> 	<p>Difficulty understanding loud speech. Total communication may be required</p> 
In-the-Ear (ITE) type Ready-made (Digital) For Right Ear 	HC-A1			
In-the-Ear (ITE) type Ready-made (Digital) For Left Ear 	HC-A1			
In-the-Ear (ITE) type Ready-made 	HM-04			
Behind-the-Ear (BTE) type 	HB-23P-24			
	HB-37S-37T			
			HB-79P	
Body-worn type 	HA-20DX			
			HA-27DX	