

Sphygmotron™

JAPAN

MERCURY-FREE SPHYGMOMANOMETER

LCD DISPLAY



Model No: ZW-302

Specifications:

- Measure unit:mmHg
- Minimal scale:LCD column:2mmHg
- Numerical display:1mmHg
- Measure method:stethoscope
- Measure scope:0-300mmHg
- Accuracy: +/-3mmHg
- Pulse rate:30-200m,+/-5%
- Pressurization:manual by bulb
- Depressurization:Manual by air release valve
- Power supply:4.5V,AA

H.SA. Device Reg. No: DE0018392

Authorised Importer & Distributor:

ZOWA
(Since 1977)

Zohra Wahab Trading Company Pte Ltd , Singapore

TEL: (65)-6743 6050 / 6747 0329

URL: www.zowamedical.com.sg