

XJ-2002DS Arm Type

DESCRIPTION

- Automatic measurement of systolic, diastolic and pulse
- Intelligent inflate, reducing the pressure on the arm
- Big LCD screen easy to read
- 80 measurement results with time & date
- Indication of WHO blood pressure standard
- Automatic power off in 3 minutes
- Modern and compact design



KD-202A (XJ-2002DS)	
FULL-AUTOMATIC ARM BLOOD PRESSURE MONITOR	
KEY SPECIFICATIONS	
Measurement method:	Oscillometric
Display type:	Digital LCD 45mm [L] x 58mm [W]
Measuring range:	Pressure 0-300 mmHg(0-40kPa)
	Pulse 40-199 p/min
Measuring accuracy:	Pressure $\leq \pm 3$ mmHg (± 0.4 kPa)
	Pulse $\leq \pm 5\%$
Operating environment:	Temperature 5°C-40°C
	Relative humidity $\leq 80\%$
Reservation environment:	Temperature -40°C-55°C
	Relative humidity $\leq 93\%$
Power source:	DC6V 3W, 4 x AA batteries
Dimension:	138mm [L] x 112mm [W] x 60mm [H]
Weight:	About 280g
PACKING DETAILS:	
Color box size:	162 mm [L] x 105mm [W] x 128mm [H]
Carton size:	560mm [L] x 355mm [W] x 300mm [H]
Packing:	20pcs/ctn

Blood Pressure Classification Chart

The WHO BLOOD PRESSURE CLASSIFICATION CHART is a general guideline for individual's blood pressure range. However, each individual's blood pressure varies with ages, circumstances... It is important to understand that this chart is for reference instead of clinic uses, please consult with your physician regularly.

