

Sphygmotron™

JAPAN

Arm-type Fully Automatic
Blood Pressure Monitor

ZW-1359b



Japan Algorithm Technology

Clinically Proven Accuracy

Detects Irregular Heart Beat

60 x 2 Memory Groups

Large Cuff (22 - 42 cm)

Talking Mode

Easy Carrying Pouch

**FREE AC
ADAPTOR!**

**LARGE LCD
WITH
BACKLIGHT**



CE 0197

Importer & Distributor:

ZOWA

Tel: (65)-6743 6050 / 6747 0329

URL: <http://www.zowamedical.com.sg>

Sphygmotron™

JAPAN

Arm-type Fully Automatic Blood Pressure Monitor

ZW-1359b



LARGE LCD WITH BACKLIGHT

Japan Algorithm Technology

Clinically Proven Accuracy

Detects Irregular Heart Beat

60 x 2 Memory Groups

Large Cuff (22 - 42 cm)

Talking Mode

Easy Carrying Pouch



Product Specification

Measurement Method	Oscillometric Method	
Measurement Range	Systolic pressure	60mmHg~260mmHg
	Diastolic pressure	30mmHg~200mmHg
	Pulse	30 to 180 Beats/Minute
Measurement Accuracy	Pressure	±3 mmHg
	Pulse	±5%
Pressurization	Automatic Pressurization	
Memory	120 Memories in Two Groups with Date and Time	
Function	Irregular Heartbeat Detection	
	WHO Classification Indicator	
	Last 3 Results Average	
	Low Battery Detection	
	Automatic Power-Off	
	Voice	
	backlight	
Power Source	4 x AAA batteries	
Unit Weight	Approx. 480g (16.93oz.) (Excluding Battery)	
Unit Dimensions	Approx. 131.2×101.8×44.1mm (5.17" x 4.01" x 1.74")(L x W x H)	
Operating Environment	Temperature	10°C ~ 40°C (50°F~104°F)
	Humidity	15% ~ 93% RH
	Pressure	700 hPa-1060 hPa
Storage Environment	Temperature	-25°C~70°C (-13°F~158°F)
	Humidity	≤ 93% RH

Importer & Distributor:

ZOWA

Tel: (65)-6743 6050 / 6747 0329

URL: <http://www.zowamedical.com.sg>

CE 0197