



WALKERS COMPANION®

AXILLARY CRUTCHES



SPECIFICATIONS:

- PUSH BUTTON FOR HEIGHT ADJUSTMENT - 9 HOLES
- ANODIZED ALUMINUM
- PADDED AXILLARY SUPPORT
- WORKING LOAD - 100 - 150KGS
- PADDED SUPPORT & HANDLE GRIP - IN GREY COLOUR

● SIZES -	<u>LENGTH</u>	<u>WIDTH</u>
SMALL /ADULT :	92 - 112 cm	20 cm (4ft 6" > 5ft 2")
MEDIUM /ADULT :	102 - 122 cm	20 cm (5ft 2" > 5ft 10")
LARGE / ADULT/TALL :	131 - 151 cm	20 cm (5ft 10" > 6ft 6")

Imported & Distributed by:

Zohra Wahab Trading Co Pte Ltd
629 Aljunied Road, #08-02
Cititech Industrial Road, #08-02
Tel: 65-6743 6050 / 65-67471806